

WormScanner Installation Notes

WormScanner is part of [QuantWorm](#) (a software package for automatic phenotyping in *C. elegans*).

Installation Instructions

It was verified to work using Java version 1.6 (or higher)

1. Download the zip file of WormScanner and extract it. A folder named WormScanner will contain the software.
2. WormScanner uses video and stage components from WormTracker 2. Due to their license restrictions, please download their software directly from their [website](#).
3. Install WormTracker by following their installation instructions.
4. After installing WormTracker, find the main 'jar' file of WormTracker, likely with file name: WormTracker.jar
5. Copy such jar file to the 'lib' folder of WormScanner
6. Compile WormScanner
 - We utilize ANT to compile WormScanner; there is a build.xml already inside the WormScanner folder.
In case that you get an error about missing/unknown CameraInfo, then replace the two occurrences of CameraInfo.CMU1394 by "CMU 1394 Digital Camera Device"
Note that the 1394-camera needs to be setup to be recognized via Java's JMF, such as following these steps: [camera_1394](#)
 - Otherwise, compile WormScanner by including into the classpath each of the jar files inside the 'lib' folder.
7. Run the program via clicking on the file: runit.bat

[Zhong's lab](#).

Contact: Boanerges Aleman-Meza (mag AT wormloco.org).

Last update on 12/18/2016